

® mod. dép. CE

SWISS MADE



# LEGHOLDERRX.CH

THE „MUST“ FOR SURGERY OF THE KNEE JOINT

LegholderRX® LEG POSITIONERS FOR KNEE SURGERY  
FLYER 2021

SWISS MADE  SINCE 2000



**novus med**

SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

## GENERAL INFORMATION



6 different stable positions. In combination with a lateral post at the middle of the thigh, one assistant becomes superfluous. X-ray transparent. Get fixed under the sterile towel.

## TECHNICAL DATA

Dimensions	490 x 240 x 200mm
Weight	3 kg, including 2 robust (long or short) velcro belts
X-ray	transparent
Material	Low maintenance, robust, non-allergenic plastic
Cleaning	No need to sterilize, inert against common cleaning agents



T: +44 (0) 1926 895089  
E: [customerservice@novus-med.com](mailto:customerservice@novus-med.com)  
W: [novus-med.com](http://novus-med.com)



**novus med**

SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

© 2024 novus med. All rights reserved. Company Registration: 10242963



## IMPROVED STRAPS

The straps have a practical locking bar and an integrated drawing for correct use.

Available in 2 lengths:

- short straps (ideally used when you can fasten the LegholderRX around one leg section)
- long straps can be used over one leg section as well as over both leg sections.



SWISS MADE  SINCE 2000

- LegholderRX®
- mod. dép.
- CE (Swissmedic CH-201505-0035)
- FDA registered
- GLN: 7649994277000

T: +44 (0) 1926 895089  
E: [customerservice@novus-med.com](mailto:customerservice@novus-med.com)  
W: [novus-med.com](http://novus-med.com)



**novus med**

SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

© 2024 novus med. All rights reserved. Company Registration: 10242963