



Imperative for
easy access and
safe procedures

novus med
SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

bioradmedisysTM
science for people





Features:

Ergonomic hub design


for easy one hand use

Tapered dilator tip (6 Fr)

for easier intramural access

Iridium / Platinum marker band at the sheath proximal end

for enhanced radiopacity / visibility under fluoroscopy, facilitating precise positioning



Kink resistant sheath, reinforced with SS

for better stability and strength while navigating through tortuous anatomy

Also available with

Hydrophilic coating

for ease of insertion

UROX™ Advantages:

Reduces the chances of trauma to the ureter during repetitive instrument exchanges.

The continuous working channel also protects delicate instruments and flexible endoscopes from damage, avoiding costly repairs.

Ordering Information	Length / Size	9.5 / 11.5 Fr	10.7 / 12.7 Fr	12 / 14 Fr
Ureteral Access Sheath (with Hydrophilic coating)	35 cm	BUUAS0309535	BUUAS0310735	BUUAS0312035
	45 cm	BUUAS0309545	BUUAS0310745	BUUAS0312045
	55 cm	BUUAS0309555	BUUAS0310755	BUUAS0312055

novus med
SURGICAL INNOVATION. KNOWLEDGE & RELIABILITY

novus-med.com
Novus Med Ltd

Registered Office
1110 Elliott Court Coventry Business Park, Herald Avenue, Coventry, West Midlands, United Kingdom, CV5 6UB. T: +44 (0) 1926 895 089

bioradmedisys™
science for people

biorad-medisys.com
Biorad Medisys Pvt. Ltd.

Corporate Office
Survey No.: 48/3 & 48/7, Pashan Sus Rd, Sus Village, Taluka Mulshi, Pune - 411021. MH, India.
T: +91-020-35000200 F: +91-020-3091 2050
Bangalore Office & Manufacturing Unit
Biorad House, #660, Hulimavu, Bannerghatta Rd, Bangalore - 560076, KA, India.
T: +91-80-2648 0683 / 6565 8007 F: +91-80-2648 0684