




## RELEEZ STENT

Used for internal drainage from the ureteropelvic junction to the bladder. Highly radiopaque for better visualization and has an extremely low encrustation tendency.

\* Also available with marker pusher.

**Available in:** Open end, Closed end,  
One end multi loop, Both end multi loop




### **DOUBLE 'J' STENT - SHORT TERM**

**Sizes:** 3 - 8 Fr  
**Lengths:** 8 - 32 cm

### **BIOFLEX - MID TERM**

**Sizes:** 3 - 8 Fr  
**Lengths:** 8 - 32 cm



### **PRETOFLEX - LONG TERM**

**Sizes:** 4.5 Fr, 4.8 Fr, 5 Fr, 5.5 Fr, 6 Fr  
**Lengths:** 24 cm, 26 cm, 28 cm, 30 cm