



Single-use RhinoLaryngoscope

RH Series



Single-use



Portable



No Cross-infection



Always Ready To Go

—

For your easier
Rhinolaryngoscopy
in OPD




Specification

Model	O.D.	I.D.	Deflection	Field of View	Working Length	Direction of View
RH-M22	2.2mm	0mm	210° Up/ Down	120°	350mm	0°
RH-M32	3.2mm	1.2mm	210° Up/ Down	120°	350mm	0°
RH-M40	4.0mm	1.2mm	210° Up/ Down	120°	350mm	0°

RH Series Imaging Part

Image processing, display and recording
All-in-one



Specification

Model	VLM-03
Size	8 inches
Weight	1.05 Kg
Storage	32 GB
Display resolution	1024*768 pixels
Video output capacity	HDMI, max 1920*1080 pixels
Battery capacity	>6400 mAh
Battery working time	Maximum 4 hours
Function	Photo, video, white balance, freeze

Ordering information

Product	Item code
Single-use RhinoLaryngoscope	
RH-M22	2018B-PA005
RH-M32	2018B-PA004
RH-M40	2018B-PA008
Imaging Processor	
VLM-03	None

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.