



Single-use Choledochoscope

CH Series

Single-use



Cost-effective



Portable



No Cross-infection



Always ready to go



An ideal bedside solution of T-tube Tract Choledochoscopy



2797

An Endoscopy Solution Company

P/N: ENG-CATA-fCHS-suCH-STD-202304V1.0Rev1.0

Specification

Model	O.D.	I.D.	Deflection	Field of View	Working Length	Direction of View
CH-M32	3.2mm	1.2mm	210° Up/ Down	120°	350mm	0°
CH-M40H	4.0mm	1.2mm	210° Up/ Down	120°	350mm	0°
CH-M52H	5.1mm	2.2mm	210° Up/ Down	120°	350mm	0°

CH Series Imaging Part

Image processing, display and recording
All-in-one



Specification

Model	VLM-03
Size	8 inches
Weight	1.05 Kg
Storage	32 GB
Display resolution	1024*768 pixels
Video output capacity	HDMI, max 1920*1080 pixels
Battery capacity	>6400 mAh
Battery working time	Maximum 4 hours
Function	Photo, video, white balance, freeze

Ordering information

Product	Item code
Single-use Choledochoscope	
CH-M32	2019B-PA004
CH-M40H	2019B-PA005
CH-M52H	2019B-PA006
Imaging Processor	
VLM-03	2021B-PA003

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.