



Single-use RhinoLaryngoscope

RH Series



Single-use



Portable



No Cross-infection



Always Ready To Go



—

For your easier
Rhinolaryngoscopy
in OPD



An Endoscopy Solution Company


P/N: ENG-EU-CATA-frHS-suRH-INTENT-202306V1.0

Specification

Model	O.D.	I.D.	Deflection	Field of View	Working Length	Direction of View
RH-M22	2.2mm	0mm	210° Up/ Down	120°	350mm	0°
RH-M32	3.1mm	1.2mm	210° Up/ Down	120°	350mm	0°
RH-M40	4.0mm	2.0mm	210° Up/ Down	120°	350mm	0°
RH-M50	5.0mm	2.8mm	210° Up/ Down	120°	350mm	0°
RH-M52	5.0mm	2.6mm	210° Up/ Down	120°	350mm	0°
RH-M58	5.8mm	3.4mm	210° Up/ Down	120°	350mm	0°

RH Series Imaging Part

Image processing, display and recording
All-in-one



Specification	
Model	VLM-03
Size	8 inches
Weight	1.0 Kg
Storage	32 GB
Display resolution	1024*768 pixels
Video output capacity	HDMI, 1920*1080 pixels
Battery capacity	≥6400 mAh
Battery working time	Up to 4 hours

Ordering information

Product	Item code
Single-use RhinoLaryngoscope	
RH-M22	2018B-PA005
RH-M32	2018B-PA004
RH-M40	2018B-PA003
RH-M50	2018B-PA001
RH-M52	2018B-PA002
RH-M58	2018B-PA006
Image processor	
VLM-03	2021B-PA003

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

HugeMed

Shenzhen HugeMed Medical Technical Development Co., LTD.

6F, Block A, Tempus Building, 1st Qingshuihe Road, Luohu District, Shenzhen City, Guangdong, China.

+86-755-22275833

+86-755-22275833

info@hugemed.net

For more information
please visit our website:
www.hugemed.net