



Single Use Solutions for Outpatient Flexible Cystoscopy

We offer a complete flexible cystoscopy solution, with competitive pricing on the single use cystoscope, monitor system and instrumentation for stent removal, biopsy and Botox injection.



HD Flexible Scope: A 15.5Fr High-Definition scope with 6.6Fr working channel, integrated suction and 210° up/down deflection is ideal for the outpatient department when combined with our 15" monitor. A portable 8" monitor is also available with a 12Fr scope for difficult catheterisations.



S-Curve Dilators: Designed for atraumatic dilation of the male urethra, these hydrophilic dilators come in a set of seven dilators between 8-20 Fr. Individually packed 22Fr, 24Fr & 26Fr dilators are available separately.



Grasping & Biopsy Forceps: Rat Tooth-Alligator Forceps are ideal for stent removal and biopsy cup forceps for bladder biopsy. The 1.8mm diameter ensures a good irrigation flow during use and the 600 mm working length is optimised for flexible cystoscopes.



Deployable Injection Needle: The 4mm, 25G diamond cut needle is suitable for injections into the bladder wall for Overactive Bladder and Neurogenic Detrusor Overactivity. The needle is housed within the sheath for safe deployment.

Product Code	Size
CY-M52H	15Fr HD flexible cystoscope, 6.6Fr working channel. Box of 10
BUSD0382040H	S-Curve urethral dilator set 8-20Fr, 40cm, Hydrophilic coating
EGF-C-1806-N	Rat Tooth-Alligator Forceps 1.8mm x 600 mm. Box of 25
EBF-A-1806	Oval Cup Biopsy Forceps 1.8mm x 600 mm. Box of 25
ATE-ZSZ-18x1200x25x4	25G deployable injection needle, 1.8mm x 1200mm. Box of 10

