




HUGEMed

novus med

SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

7.5Fr & 9.0Fr Single-use Flexible Ureterorenoscopes



Brilliance: Excellent visualisation due to the latest CMOS sensor with 160,000-pixel resolution and dual distal LED illumination.


Integration: The HUV-01 control unit integrates and is viewed through your existing camera system. Images and video can be captured from the scope and saved to an external USB storage device.



Functionality: Both the 7.5Fr and 9Fr models offer 285° up-down deflection, 1:1 torque stability and a standard 3.6Fr working channel. Both European and US deflection versions are available.

Innovation: Your choice of 3 o'clock or 9 o'clock instrument channel positions is available, optimising laser fibre access dependent on stone location:



	Right Kidney	Stone Position	Left Kidney
	3 o'clock	Lower Pole	9 o'clock
	3 o'clock	Posterior Calyx	9 o'clock
	9 o'clock	Anterior Calyx	3 o'clock
			

Deflection	Channel	Product Code 7.5Fr	Product Code 9.0Fr
Standard / US	9 o'clock	HU30SSL	HU30SL
	3 o'clock	HU30SSR	HU30SR
Reverse / Europe	9 o'clock	HU30SRL	HU30RL
	3 o'clock	HU30SRR	HU30RR



novus med
SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

LG1.SUFU.V1

+44 (0) 1926 895089
customerservice@novus-med.com

