



Mediflex®

novus med
SURGICAL INNOVATION, KNOWLEDGE & RELIABILITY

Mediflex® Cook®-Style Nathanson Liver Retractor

Mediflex® Nathanson Liver Retractors provide atraumatic retraction for laparoscopic and robotic procedures. Their secure hexagonal end prevents slippage. Mediflex® offers the largest selection, suitable for all patients, including bariatrics and paediatrics.




Ease of Insertion: an extended tip length is available for ease of insertion through thicker abdominal walls.



Low-Profile: Robotic version has low-profile design with shorter shaft to reduce the risk of collision with robot arms.



Secure: Direct attachment to Mediflex® StrongArm® and FlexArm™ systems prevents liver slippage and axial rotation.



Product Code	Description	Length (mm)					
		A	B	C	D	E	F
69735-C	Small 6.5mm Cook®-Style Retractor	22	14	50	44	89	140
69735-C-XL	Small 6.5mm Cook®-Style Retractor, Long Tip	40	14	50	44	89	140
69736-C	Medium 6.5mm Cook®-Style Retractor	25	20	51	52	127	196
69736-C-R9	Medium 6.5mm Nathanson Retractor, Cook®-Style, Robotic Shaft 9cm	25	20	51	52	127	51
69736-C-XL-R9	Medium 6.5mm Cook®-Style Retractor, Long Tip, Robotic Shaft 9cm	44	20	51	52	127	51
69736-C-XL	Medium 6.5mm Cook®-Style Retractor, Long Tip.	44	20	51	52	127	196
69737-C	Large 6.5mm Cook®-Style Retractor	38	27	76	70	178	272
69737-C-XL	Large 6.5mm Cook®-Style Retractor, Long Tip	57	27	76	70	178	272