


The New CS Universal Spine Extension

The CS Universal Spine Extension can be fitted to a **wide range of standard operating tables** offering a patient position and imaging capabilities **comparable to that of a specialist dual-rail spine table**.


- + Attaches to **ANY** standard operating table
- + **Small Storage Footprint**
- + **360° C-Arm/O-Arm access**
- + **275kg/600lbs. patient weight capacity**
- + **Hinge mechanism to allow flexion and extension**
- + **Compatible with Mayfield/Doro skull clamp**



**CS Universal Spine
Extension - EU**

Code: CSM2958
Qty: 1 Set

Whats Included:



CS Universal Spine Extension (Set)

Code: CSM2958
Qty: 1 Set

Includes:

CSM2957 - CS Universal Spine Extension - Frame
CSM2955 - CS Universal Spine Extension - Table Adaptor (EU)
CSM2956 - CS Universal Spine Extension - Powered Leg



CS Prone Plus Head Support System

Code: CSM2565
Qty: 1 System
(Helmet & Mirror)



CS Universal Arm Support - Moulded

Code: CSM2647
Qty: 1 Unit



CS Easy Shin Support

Code: CSM3020
Qty: 1 Unit



CS Knee Rings (Pair)

Code: CSM2821
Qty: 1 Pair



CS Table Rail Adaptor - EU

Code: CSM2739-1
Qty: 1 Unit



CS Prone Plus Face Cushion

Code: CSM2512
Qty: 10 Units



CS Head Tray Attachment - Lowered

Code: CSM2751
Qty: 1 Unit



CS Chest Support

Code: CSM2709
Qty: 1 Unit



CS Hip Supports (pair)

Code: CSM2717
Qty: 1 Pair



CS Thigh Supports (pair)

Code: CSM2718
Qty: 1 Pair



CS Pad 2"

Code: CSM2568-2
Qty: 1 Unit



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Units

Optional Accessories:



CS Din Rail Clamp

Code: CSM2763
Qty: 1 Unit



CS Supine Top-Base Only

Code: CSM2722-7
Qty: 1 Unit



CS Supine Top - 2" Foam Pad

Code: CSM2722-2
Qty: 1 Unit



CS Comfort Covers - 5 Piece

Code: CSM2600
Qty: 6 Units



CS Carbon & Universal Extension Cart

Code: CSM2954
Qty: 1 Unit

Contact Us For More Information

Telephone: +44(0) 1926 895 089

Email: customerservice@novus-med.com

Web: novus-med.com

