

DISPOSABLE BOUGIE DILATOR SETS

Intended use

For dilation and treatment of the oesophagus and cardia

Key features

- **Medical-grade material for comfort and flexibility**
- **Secure anti-shedding design for patient safety**
- **Smooth tip for easy, atraumatic insertion**
- **Clear X-ray visibility with radiopaque markers**
- **Printed scales for accurate depth control**



Model	Dilator diameter (mm)	Working Length (mm)	Guide wire
BDS-05	5	850	0.035"
BDS-06	6	850	0.035"
BDS-07	7	850	0.035"
BDS-08	8	850	0.035"
BDS-09	9	850	0.035"
BDS-10	10	850	0.035"
BDS-11	11	850	0.035"
BDS-12	12	850	0.035"
BDS-13	13	850	0.035"
BDS-14	14	850	0.035"
BDS-15	15	850	0.035"
BDS-16	16	850	0.035"
BDS-17	17	850	0.035"
BDS-18	18	850	0.035"
BDS-19	19	850	0.035"
BDS-20	20	850	0.035"