


HugeMed Single-use Flexible Cystoscopes



Portable: Ideal for difficult catheterisations, the 15.6" HD display monitor can be easily moved between departments or for a more permanent solution, our HUV-02 control unit can be integrated into an existing camera system in theatres or outpatients.



TULA Compatible: Full compatibility with both 980 nm and 1470 nm wavelengths, avoiding the interference seen with some other manufacturers' scopes.



Irrigation Control: A Tuohy-Borst with side-arm adapter and stopcock is included with every scope. This is cost saving and enables controlled irrigation from the sterile field.



Suction as Standard: The suction function and large 2.8mm channel (CY-M50) helps clear the visual field by removing fluid and debris providing better visualisation of the area being examined during cystoscopy.



Choice of Sizes: Available in the standard 15Fr version with a \varnothing 2.8mm / 8.4Fr channel or a slimmer 12Fr scope for tighter anatomies, patient comfort or as a flexible mini-nephroscope with a \varnothing 2.0mm / 6.0Fr channel.

Product Code	Description	NPC
CY-M52H	15Fr HD flexible cystoscope, \varnothing 2.2mm working channel, box of 20	FBE88300
CY-M50	15Fr flexible cystoscope, \varnothing 2.8mm working channel, box of 20	FBE88301
CY-M40	12Fr flexible cysto-nephroscope, \varnothing 2.0mm / 6.0Fr working channel, box of 10	FBE86171

