



BARI calibration bougies

AVAILABLE IN 34, 36, 38, 40 AND 46 FRENCH

- Calibration tube for administration of meth-blue to check the integrity of the staple line is located at the distal end of the catheter where it is needed. The tube for administering meth-blue is 2mm in diameter. This means less wastage and administration at desired site and avoids any unintentional leakage.
- Dual Lumen, made from medical grade silicone with atraumatic design for easy insertion and patient safety. Prevents the tube becoming weak and floppy.
- Soft conical tip to allow easy positioning in the stomach with minimal resistance.
- Distal end of the tube is soft and has holes on all sides to allow safe suction.



Straight Calibration Tube

PRODUCT CODE	SIZE	BOX SIZE
ACS2025-0034	34 French	10
ACS2025-0036	36 French	10
ACS2025-0038	38 French	10
ACS2025-0040	40 French	10
ACS2025-0046	46 French	10